



---

## INSTRUCTIONS ON HOW TO USE MERCON™ PRODUCTS IN MERCURY SPILL KIT MODEL NO. SEI268

### HOW TO USE MERCON™ PRODUCTS

1. Block off spill site (minimum 2 metre radius) from foot traffic. Open windows and shut off ventilation system. Check clothing and footwear for mercury and remove contaminated apparel.
2. Wear protective clothing, eyewear, gloves and ventilation according to applicable regulations. Spray **MERCON™ spray** into ambient air above the clean up site.
3. Use **MERCON™ vap** to cover visible mercury beads. **CAUTION:** After applying **MERCON™ spray** or **MERCON™ vap**, floor may be slippery.
4. Use **MERCURY Aspirator** to pick up mercury beads and deposit into **MERCON™ tainer**. Ensure tip is below rim.
5. **NOTE: MERCON™ tainer may be re-used, provided it is no more than 1/2 full, for up to six months. Always have a usable container available for future spill.**
6. Wipe contaminated surfaces clean with **MERCON™ vap** liquid or **MERCON™ wapes** ( kits 1 & 2 only)
7. Wipe shoes, gloves and other contaminated items with **MERCON™ wapes** or disposable towels saturated with **MERCON™ vap**.
8. Label contaminated items.
9. Determine the appropriate disposal method according to local, state and federal regulations for waste.
10. File an appropriate incident report in accordance with company or local environmental regulations. Ensure that any staff which may have come in contact with the mercury spill or which may have been exposed to mercury vapor seek medical attention as soon as possible.

**ANY QUESTIONS CALL ROSS HEALTHCARE, INC. @ 1-800-663-8303**

### MERCON™ ASPIRATOR DIRECTIONS FOR USE

The **MERCON™ Aspirator** is a simple and easy to use waste mercury pick-up tool. Effective on hard surfaces and under water. Kit allows for quick clean up of spilled mercury. When purchased in one of three Mercury Spill Kits it allows for clean up of mercury spills without any mixing or assembly.

Use of **MERCON™ Syringe Aspirator Mercury to Pick-up Mercury droplets**

1. Push loading knob all the way down.
2. Release trigger and suction up 1-2 oz. of **MERCON -VAP OR MERCON-X**. This exercise will familiarize user with action of tool.
3. Place tip of Aspirator near mercury beads and repeat steps 1 and 2.
4. Tilting tip above body of Aspirator will prevent mercury from being accidentally released.
5. Hold the syringe/aspirator above **MERCON™ tainer** being careful to ensure tip is below upper rim
6. Gently press loading knob to slowly release any mercury into container
7. Continue with spill procedure as outlined in Mercury spill kit instructions.
8. Repeat steps 1 and 2 to clean Aspirator. Wipe down exterior surfaces using one of the following **MERCON™** products: **MERCONwapes**, **MERCONvap**, **MERCON-X**